



Voss Medical Products-VM Positioning, Inc

Product #N-1248

Voss NOTCH-Nasal/Oral Intubation Tube Control Holder

Indications for use: The NOTCH is used on the forehead or the chin to hold the intubation tube in place during surgical procedures.

| <u>Specifications</u> | |
|---|-------------------|
| Color | Blue |
| Dimension_Height | 0.25 in. |
| Eco Friendly Item | No |
| HPIS Code | unknown |
| Latex Free | Yes |
| Length | 2 in. |
| Material | Polyurethane foam |
| Type of Cover | No Cover |
| UNSPSC | 42211503 |
| Width | 1 in. |
| Sterile or non-sterile? | Non-sterile |
| If non-sterile, is a sterile version available? | No |
| Contraindications | None known |